

A REVOLUTIONARY PERFORMANCE

MetalClean Plus

with **SOLIDSAFE™** technology

Mild alkaline-based highly concentrated solid cleaning agent for DI to soft water qualities

MetalClean Plus is part of the unique SolidSafe System that revolutionises the automated cleaning of surgical instruments.

Offering significant safety and ergonomic benefits, it also provides excellent cleaning performance.

- ▲ Innovative format improving safety standards in the cleaning environment
- ▲ The unique formula delivers excellent performance with reliable cleaning outcomes
- ▲ Sustainable benefits



MetalClean Plus

with **SOLIDSAFE**™ technology

Mild alkaline-based highly concentrated solid cleaning agent for DI to soft water qualities

Why Use MetalClean Plus Solid?

EXCELLENT CLEANING PERFORMANCE

A unique formula based on alkaline ingredients and complexing agents, MetalClean Plus is low foaming and ensures excellent cleaning performance with reliable cleaning outcomes. The formula is designed to be used in DI water to soft water.

The good material compatibility gives maximum protection to your valuable instrumentation during reprocessing.

REVOLUTIONARY FORMAT

The highly concentrated solid chemistry, in an easy to handle 4 kg capsule, provides the equivalent of 80L of liquid which means:

- ▲ Less weight
- ▲ Less space for storage
- ▲ Less waste

INNOVATIVE SYSTEM

Using the solid in a closed dispenser means that there will be no risk of leaking and spilling canisters and drums. The revolutionary SolidSafe System does away with the messy changing of dip tubes or suction lances.

The coded capsules for the appropriate dispenser ensures that there is no switched chemistry.

ENVIRONMENTALLY RESPONSIBLE

MetalClean Plus Solid is free from phosphates, EDTA and NTA. The packaging can be disposed of in the regular waste which gives a great reduction in disposal fees, with no chemical waste to harm the environment.



MetalClean Plus 4 kg Capsule



SolidSafe Dispenser

SIMPLE. INNOVATIVE. SAFE.

APPLICATION AREA

MetalClean Plus is suitable for instruments and utensils made of stainless steel, instrument steel, optics, conventional plastics and materials of anesthesia equipment. MetalClean Plus has been extensively tested on anodized aluminium, however it is recommended to test your specific anodized instruments prior to the use for confirmation of suitability.

INSTRUCTIONS FOR USE

MetalClean is suitable for instruments and utensils made of stainless steel, instrument steel, optics, conventional plastics and materials of anaesthesia equipment.

Recommendations:

3-10ml/L of the solution produced by the SolidSafe Dispenser is automatically drawn via the alkaline resistant dosing pumps of the washer disinfectant. Select program according to washer manufacturer's instructions.

INGREDIENTS

MetalClean Plus: Sodium hydroxide, complexing agents and auxiliary substances.

PACKAGING AND ORDERING INFORMATION

PRODUCT	OUTER SIZE	ORDER CODE
MetalClean Plus Solid Capsule	2 x 4kg	3071060

STORAGE

Capsules should be stored in a cool, dry location out of direct sunlight. To minimise lifting, store on a shelf at waist height.

ECOLAB GULF FZE

Tel. : +971 4 801 4000

Fax : +971 4 884 1782

P.O.Box : 262015 JAFZA, DUBAI, UAE

E-mail : ecolab.healthcaremea@ecolab.com

www.ecolab.com



ECOLAB®